

Persona IQ[®]

THE SMART KNEE[™]






Introducing **Persona IQ Knee Implant** The Smart Knee

The size, shape, and fit of your knee implants can make a big difference in your satisfaction after total knee replacement. That's why the Persona Knee was designed to give your surgeon more options to better match your bone's anatomy.

We are proud to introduce the first and only smart knee implant, Persona IQ. Through the latest implant technology, you'll have confidence that you're connected to your healthcare team, every step of the way.



SCAN HERE
to see our
patient video



At Zimmer Biomet,
we're never done innovating and
**transforming what's
possible...for you.**

WARNING - The kinematic data from this device have not been demonstrated to have clinical benefit. It is not intended to be utilized for clinical decision-making, and no data have been evaluated by FDA regarding clinical benefits.

How does it work?

Your surgeon will start by selecting Persona Knee implants that best match your knees' anatomy. Then, on the bottom portion of the implant, a small "smart" stem extension is attached. This stem contains sensors that capture information like your range of motion, stride length, walking speed, and step count. You, your surgeon, and medical care team will be able to view your data via a secure online portal.



Why collect data?

The information collected by Persona IQ allows your physician to remotely view your knee data throughout your recovery. Through the latest implant technology, you'll have confidence that you're connected to your healthcare team, every step of the way.



Sounds great, right? Well, what if we told you that Persona IQ is easy to use too!

Persona IQ fits seamlessly into your life. You'll simply set it up one time, following instructions that will arrive in your welcome kit. Then, it automatically collects data overnight. This data is sent directly to your surgeon.

Persona IQ collects data on:



The range of motion
in your knee



Your step count



Your stride length



The distance you walk



Your walking speed

Persona IQ does NOT collect data on:



Your location





What happens next?

If you're interested in Persona IQ, The Smart Knee, talk to your doctor as this product is not yet available in all areas. Once you and your surgeon have decided Persona IQ is right for you, they will enroll you and guide you on how to get started.

Your welcome kit will include:

- Quick start guides
- User manual
- Base station with USB cable and wall adapter
- Patient brochure



We care about your privacy



Privacy is a fundamental human right and Zimmer Biomet maintains robust practices to help ensure that right is protected. Zimmer Biomet will only use customer and patient data to provide the services we have agreed upon, and for purposes that are compatible with providing those services. Zimmer Biomet protects the data that you entrust us with through strong safeguards.



To learn more about Persona IQ,
email SmartKneeSupport@zimmerbiomet.com
or call 844-799-8208.

For product information,
talk to your surgeon and visit www.zimmerbiomet.com.

Legal Manufacturer:

Zimmer, Inc.
1800 West Center St.
Warsaw, Indiana 46581-0587
USA
zimmerbiomet.com

Legal Manufacturer:

Canary Medical USA LLC
2710 Loker Ave, West, Suite 350
Carlsbad, CA 92010 USA

Disclaimers

This brochure is intended to provide an overview of knee replacement surgery and should be reviewed with your doctor. It does not include all of the information needed to determine eligibility for knee replacement or for the proper use and care of artificial knee replacements. Please consult your surgeon for more information. Information may also be obtained by calling the toll-free number or visiting the web site. The toll-free number also can be used to obtain product contraindications, warnings, precautions, and possible adverse effects. Individual results may vary. Your results will depend on your personal circumstances. How long a knee replacement will last varies from patient to patient. It depends on many factors, such as the patient's physical condition, activity level, and body weight and the surgical technique. Replacement joints are not as strong or durable as a natural, healthy joint, and there is no guarantee that an artificial joint will last the rest of a patient's life. All knee replacements may need to be replaced at some point. All content herein is protected by copyright, trademarks and other intellectual property rights, as applicable, owned by or licensed to Zimmer Biomet or its affiliates unless otherwise indicated, and must not be redistributed, duplicated or disclosed, in whole or in part, without the express written consent of Zimmer Biomet.

©2021 Zimmer Biomet.

