ACD-A

Anticoagulant Citrate Dextrose Solution,
Solution A, U.S.P., ACD-A, PN 9350, 30 mL
NDC 23731-6051-3

Description:

Anticoagulant Citrate Dextrose Solution, Solution A, U.S.P., (ACD-A), is a sterile, non-pyrogenic solution of citric acid, sodium citrate, and dextrose, in water for injection. Each 10 mL of solution contains:

Citric Acid, anhydrous, USP	0.073 g
Sodium Citrate, dihydrate, USP	0.220 g
Dextrose, monohydrate, USP	0.245 g
Water for Injection, USP	q.s.

pH: 4.5 – 5.5

⊗ Single patient use only, on a single occasion.

Clinical Pharmacology:

ACD-A acts as an anticoagulant by the action of the citrate ion chelating free ionized calcium, thus making calcium unavailable to the coagulation system.

FL9350 01/20

Indications and Usage:

Anticoagulant Citrate Dextrose Solution, Solution A, U.S.P. (ACD-A), is intended for use as an anticoagulant in the extracorporeal blood processing with Autologous PRP Systems in production of platelet rich plasma (PRP). Refer to the manufacturer's Operator's Manual of the Autologous PRP System for the Directions for Use.

Contraindications:

NOT FOR DIRECT INTRAVENOUS INFUSION

Warnings:

NOT FOR DIRECT INTRAVENOUS INFUSION

Precautions:

General

- Aseptic technique must be maintained at all times.
- ACD-A is a clear and colorless solution. If the product shows any cloudiness or turbidity, the product should be discarded.
- The closure system provides a biological barrier and should be intact – discard product if system is compromised.

Information for Patients

None.

FL9350 01/20

Carcinogenesis, mutagenesis, impairment of fertility

Long-term studies in animals have not been performed to evaluate the carcinogenic potential of ACD-A.

Pregnancy

Long-term studies in animals have not been performed to evaluate the effects of ACD-A on pregnant women.

Pediatric Use

The safety and effectiveness of ACD-A in children have not been established.

Adverse Reactions:

Not applicable. This product is used as an anticoagulant in the extracorporeal blood processing with Autologous PRP Systems in production of platelet rich plasma (PRP).

Overdosage:

Not applicable. This product is used as an anticoagulant in the extracorporeal blood processing with Autologous PRP systems in production of platelet rich plasma (PRP).

Dosage and Administration:

Refer to the manufacturer's Operator's Manual of the Autologous PRP System for the Directions for Use.

Rx Only:

Caution: Federal Law (US) restricts this device to sale by or on the order of a Physician.

FL9350 01/20

How Supplied:

Anticoagulant Citrate Dextrose Solution, Solution A, U.S.P., (ACD-A)

REF	SIZE	CASE
PN 9350-10	30 mL	10 Vials/Case
PN 9350-25	30 mL	25 Vials/Case

It is recommended that the product be stored at ambient room temperature, 24 °C (75 °F); however, the product can be stored between 15 °C (59 °F) and 30 °C (86 °F). Protect from freezing and exposure to excessive heat should be minimized.

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Citra Labs, LLC
Global Supply Chain Center B.V.