

Omni™ Surgical Suite System Release Notes

Omni Suite v1.0.13 – 06/30/2021	
1.0.13	Release notes - Omni Zimmer - Version OmniSuite_Bugfix_1.0.13
1.0.13	Fixed user interface issues and addressed verbiage issues

Known Issues

- With manual workflow, if the workflow is advanced from anesthesia start to surgery start by mistake and then reset to anesthesia start, the anesthesia timer will not reset. If it is advanced from surgery start to surgery stop by mistake and then reset to surgery start the surgery timer will not be reset.

Omni Suite v1.0.16 - 08/08/2022	
1.0.14	Release notes - Omni Zimmer - Version OmniSuite_Hotfix_1.0.14
1.0.14	Updated application versions to latest .NET Framework

Omni Suite v1.0.17 - 12/12/2022	
1.0.15	Release notes - Omni Zimmer - Version OmniSuite_Hotfix_1.0.15
1.0.15	Incorporated improvements for cybersecurity
1.0.6	Release notes - Omni Zimmer - Version OmniSuite_AI_Hotfix_1.0.6
1.0.6	Fixed log file deletion issue

Omni Suite v1.0.18 - 12/22/2022	
1.0.7	Release notes - Omni Zimmer - Version OmniSuite_AI_Hotfix_1.0.7
1.0.7	Improved Exception handling