



Comprehensive® Total Shoulder System

BIOMET®

Comprehensive® Total Shoulder System

Anatomic engineering, surgeon versatility and ease of use

The Comprehensive® Total Shoulder system is an evolutionary design based on the successful clinical heritage of the Bio-Modular® Shoulder System. Utilizing innovative technologies, the system is able to allow for surgeon efficiency through intraoperative flexibility, unmatched humeral stem options and infinite humeral head offset possibilities.

Proximal PPS® Porous Plasma spray allows for enhanced biologic fixation

Fully convertible platform stem portfolio, featuring Micro, Mini, Standard and Revision lengths



Humeral Stems

- 45°/135° anatomic neck-shaft angle
- PPS® Porous Plasma Spray coating proximally to enhance biologic fixation
- Reverse Morse taper allows for increased access to the glenoid and conversion to the Comprehensive® Reverse Shoulder
- Available in Micro, Mini, Standard and Revision lengths
- Uncemented or cemented use



Versa-Dial® variable offset humeral head allows for proper anatomic reconstruction

Hybrid® Modular Glenoid allows for three variations in glenoid fixation

Regenerex® central peg allows for long-term biologic fixation

Versa-Dial® Humeral Heads

- Infinite offset options between 0.5 and 4.5 mm
- Allows for accurate reproduction of natural anatomy
- 17 head sizes, ranging from 38 to 58 mm in diameter and 18–27 mm in height



Hybrid® Glenoids

- Non-conforming diameter of curvature provides the ability to use any size head with any size glenoid
- Three outer pegs for cemented fixation
- Two central peg options
 1. Regenerex® Porous Titanium Construct for biologic fixation
 2. Polyethylene compression fit for cemented use
- Small, Medium and Large glenoid base sizes

All trademarks herein are the property of Biomet, Inc. or its subsidiaries unless otherwise indicated.

Comprehensive[®], PPS[®], Versa-Dial[®], Bio-Modular[®] and Regenerex[™] are trademarks of Biomet Manufacturing Corp.

This material is intended for the Biomet Sales force and physicians only and is NOT intended for patient distribution. It is not to be redistributed, duplicated or disclosed without the express written consent of Biomet.

For product information, including indications, contraindications, warnings, precautions and potential adverse effects, see the package insert and Biomet's website.



One Surgeon. One Patient.[®]

P.O. Box 587, Warsaw, IN 46581-0587 • 800.348.9500 x 1501
©2011 Biomet Orthopedics • biomet.com
Form No. BOI0048.1 • REV121511