

ALL-SUTURE ANCHORS ALL-INCLUSIVE P O R T F O L I O



2009

Biomet releases the 1.4 mm JuggerKnot® All-Suture Anchor, the first all-suture anchor on the market.



2011

JuggerKnot 1.4mm Short




2013

JuggerKnot launches in the Asia Pacific market. JuggerKnot 1.45mm Long



JuggerKnot launches in the EMEA market.

2010



JuggerKnot 2.9mm and JuggerKnot 1.0mm Mini

2012

THE LEADER IN ALL-SUTU

In 2009 Zimmer Biomet became the first company to introduce an all-suture anchor to the market with the launch of the JuggerKnot, revolutionizing soft tissue repair. Driven by a commitment to patient rehabilitation, Zimmer Biomet has expanded its portfolio into an array of devices to meet surgeon preferences.

2014

Cayenne Medical introduces the SureLock® anchor; Marking a new, mechanically deployed, all-suture option.



2019



Zimmer Biomet acquires Cayenne Medical, furthering the Zimmer Biomet Sports Medicine All-Suture Anchor product offerings.

2016



JuggerKnotless and JuggerKnot 1.45mm Rigid

2014



JuggerKnot with BroadBand Tape

2021


RE ANCHOR INNOVATION

JUGGERKNOT 2.9 MM ALL-SUTURE ANCHORS

The JuggerKnot 2.9 mm Soft Anchor has multiple offerings and designs to accommodate surgeon preference in arthroscopic and mini-open procedures.

LOADED WITH MAXBRAID SUTURE


JuggerKnot 2.9mm Anchors are loaded with MaxBraid Suture. The JuggerKnot 2.9 mm Anchors with MaxBraid Sutures are offered double loaded with and without needles for arthroscopic and mini-open repairs.



LOADED WITH BROADBAND TAPE

The JuggerKnot 2.9 mm Anchors loaded with BroadBand Tape are offered in a variety of designs: double and triple loaded with sliding and non-sliding sutures. Offered with and without needles for arthroscopic and mini-open repairs.

Double and triple loaded with **sliding** sutures for knotted repairs.





Double and triple loaded with **non-sliding** sutures for knotless repairs.



KNOTTED VS. KNOTLESS REPAIRS

Co-braided sutures are always sliding for arthroscopic knot tying. Solid suture colors are always non-sliding, or fixed in the anchor, designed for complete knotless repairs.



JUGGERKNOT 2.9 MM

JuggerKnot 2.9 mm with Sliding BroadBand Tape

Part Number	Description	Suture Color	Needles	QTY
CM-99129B	Double Loaded Sliding #2 (1.5 mm) BB	Wht/Blk, Blu/Wht	-	1
CM-99129BT	Double Loaded Sliding #2 (1.5 mm) BB (10 pack)	Wht/Blk, Blu/Wht	-	10
CM-99229B	Double Loaded Sliding #2 (1.5 mm) BB	Wht/Blu, Blk/Wht	-	1
CM-99229BT	Double Loaded Sliding #2 (1.5 mm) BB (10 pack)	Wht/Blu, Blk/Wht	-	10
CM-99529B	Triple Loaded Sliding #2 (1.5 mm) BB	Wht/Blk, Blu/Wht, Blk/Wht	-	1
CM-99529BT	Triple Loaded Sliding #2 (1.5 mm) BB (10 pack)	Wht/Blk, Blu/Wht, Blk/Wht	-	10
CM-99129BN	Double Loaded Sliding #2 (1.5 mm) BB with Needles	Wht/Blu, Blk/Wht	Yes	1
CM-99129BNT	Double Loaded Sliding #2 (1.5 mm) BB with Needles (10 pack)	Wht/Blu, Blk/Wht	Yes	10

JuggerKnot 2.9 mm with Non-Sliding BroadBand Tape

Part Number	Description	Suture Color	Needles	QTY
CM-99329B	Double Loaded (1) Non-Sliding and (1) Sliding #2 (1.5 mm) BB	Wht/Blu Sliding, Blk Non-Sliding	-	1
CM-99329BT	Double Loaded (1) Non-Sliding and (1) Sliding #2 (1.5 mm) BB	Wht/Blu Sliding, Blk Non-Sliding	-	10
CM-99429B	Double Loaded (1) Non-Sliding and (1) Sliding #2 (1.5 mm) BB	Wht/Blk Sliding, Blu Non-Sliding	-	1
CM-99429BT	Double Loaded (1) Non-Sliding and (1) Sliding #2 (1.5 mm) BB	Wht/Blk Sliding, Blu Non-Sliding	-	10
CM-99629B	Triple Loaded Implant (2) Non-Sliding and (1) Sliding #2 (1.5mm) BB	Wht/Blk Sliding, Blu Non-Sliding, Blk Non-Sliding	-	1
CM-99629BT	Triple Loaded Implant (2) Non-Sliding and (1) Sliding #2 (1.5mm) BB	Wht/Blk Sliding, Blu Non-Sliding, Blk Non-Sliding	-	10
CM-99229BN	Double Loaded (1) Non-Sliding and (1) Sliding #2 (1.5 mm) BB with Needles	Blu/Blk Sliding, Blk Non-Sliding	Yes	1
CM-99229BNT	Double Loaded (1) Non-Sliding and (1) Sliding #2 (1.5 mm) BB with Needles	Blu/Blk Sliding, Blk Non-Sliding	Yes	10

JuggerKnot 2.9 mm with MaxBraid Suture

Part Number	Description	Suture Color	Needles	QTY
912029	Double Loaded #2 MB	Blue/White, Blue	-	1
912529	Double Loaded #2 MB (10 pack)	Blue/White, Blue	-	10
110005093	Double Loaded #2 MB with Needles	Blue/White, Blue	Yes (Cutting)	1
110005241	Double Loaded #2 MB with Needles (10 pack)	Blue/White, Blue	Yes (Cutting)	10
110005096	Double Loaded #2 MB with Needles	Blue/White, Blue	Yes (Tapered)	1
110005242	Double Loaded #2 MB with Needles (10 pack)	Blue/White, Blue	Yes (Tapered)	10

JuggerKnot 2.9 mm Disposable Instrumentation

Part Number	Description
912091	2.9 mm Sterile Disposable Drill Bit
912057	2.9 mm Straight Instrument Kit (Guide, Obturator, Drill Bit)
110005306	2.9 mm with Needles Drill Bit Short

JuggerKnot 2.9 mm Reusable Instrumentation

Part Number	Description
912090G	2.9 mm Reusable Drill Guide
912092G	2.9 mm Reusable Obturator
912093G	2.9 mm Reusable Trocar
912064	2.9 mm Reusable In-Guide Punch
1100016951	2.9 mm Reusable Short Guide


JUGGERKNOT 1.4 AND 1.5 MM ALL-SUTURE ANCHORS

The JuggerKnot 1.4 mm All-Suture Anchor was the first of its kind, creating a revolutionary all-suture anchor market that has led to a robust family platform. For shoulder labral repair, this portfolio has expanded to include #1 and #2 MaxBraid Suture, #1 and #2 BroadBand Tape, and multiple suture color options.

**Loaded with #1
MaxBraid Suture**



**Loaded with #1
(1.0 mm wide)
BroadBand Tape**



**Loaded with #2
MaxBraid Suture**



**Loaded with #2
(1.5 mm wide)
BroadBand Tape**






#2 MAXBRAID SUTURE



#1 BROADBAND TAPE



#2 BROADBAND TAPE



JUGGERKNOT 1.4/1.5 MM

JuggerKnot 1.4 mm with MaxBraid Suture

Part Number	Description	Suture Color	QTY
CM-99114	Single Loaded #1 MB	White/Black	1
912030	Single Loaded #1 MB	Blue/White	1
912010	Single Loaded #1 MB (10 pack)	Blue/White	10
912000	Implant System (2 Implants, Straight Guide, Drill Bit)	Blue/White	10
CM-99214	Single Loaded #1 Non-Sliding MB	Black (Non-Sliding)	1

JuggerKnot 1.5 mm with MaxBraid Suture

Part Number	Description	Suture Color	QTY
CM-99115	Single Loaded #2 MB	White/Black	1
912031	Single Loaded #2 MB	Blue/White	1
912015	Single Loaded #2 MB (10 pack)	Blue/White	10
CM-99215	Single Loaded #2 Non-Sliding MB	Black (Non-Sliding)	1

JuggerKnot 1.4 mm with BroadBand Tape

Part Number	Description	Suture Color	QTY
CM-99114B	Single Loaded #1 (1.0 mm wide) BB	White/Black	1
CM-99114BT	Single Loaded #1 (1.0 mm wide) BB (10 pack)	White/Black	10

JuggerKnot 1.5 mm with BroadBand Tape

Part Number	Description	Suture Color	QTY
CM-99115B	Single Loaded #2 (1.5 mm wide) BB	Blue/Black	1
CM-99115BT	Single Loaded #2 (1.5 mm wide) BB (10 pack)	Blue/Black	10

JuggerKnot 1.4 mm Disposable Instrumentation

Part Number	Description
912036	1.4 mm Flex Drill Bit
912036R	1.4 mm Rigid Drill Bit
912136	1.4 mm Flex Drill Bit (Sterile)
912040	1.4 mm Striaight Kit (Straight Guide, Flex Drill Bit)
912040C	1.4 mm Curved Kit (Curved Guide, Flex Drill Bit)
912040P	1.4 mm Perc Kit (Perc Guide, Guide Pin, Drill Bit)
912040R	1.4 mm Straight Rigid Kit (Straight Guide, Rigid Drill Bit)
912045	1.4 mm Centering Sleeve Kit (Cntrng Slv, Flex Drill Bit)
912140C	1.4 mm Centering Sleeve System (Crvd Guide, Cntrng Slv, Flex Drill)

JuggerKnot 1.5 mm Disposable Instrumentation

Part Number	Description
912037	1.5 mm Flex Drill Bit
912037R	1.5 mm Rigid Drill Bit
912041	1.5 mm Straight Kit (Straight Guide, Rigid Drill Bit)
912041C	1.5 mm Curved Kit (Curved Guide, Flex Drill Bit)
912046	1.5 mm Centering Sleeve Kit (Cntrng Slv, Flex Drill Bit)
912141C	1.5 mm Centering Sleeve System (Crvd Guide, Cntrng Slv, Flex Drill)

JuggerKnot 1.4/1.5 mm Reusable Instrumentation

Part Number	Description
912035	1.4/1.5 mm Straight Guide
912035C	1.4/1.5 mm Curved Guide
912038	1.4/1.5 mm Rigid Trocar
912038C	1.4/1.5 mm Flex Trocar
912038P	1.4/1.5 mm Sharp Percutaneous Trocar
912039	1.4/1.5 mm Rigid Obturator
912039C	1.4/1.5 mm Flex Obturator

JUGGERKNOT 1.4 MM SHORT AND 1.45 MM RIGID ALL-SUTURE ANCHORS

The JuggerKnot 1.4 mm Short and 1.45 mm Rigid Soft Anchors were designed for guided or guideless small joint repairs depending on surgeon preference. The 1.4 mm Short All-Suture Anchors are single loaded with either a #1 MaxBraid or a #1 (1.0 mm wide) BroadBand Tape and are designed for guided insertion. The 1.45 mm Rigid All-Suture Anchors are single loaded with either #1 or #2 MaxBraid Suture or #2 (1.5 mm wide) BroadBand Tape and are designed for guideless insertion.

1.4 mm Short with MaxBraid Suture



1.4 mm Short with BroadBand Tape




1.45 MM Rigid with MaxBraid Suture



1.45 MM Rigid with BroadBand Tape





JUGGERKNOT 1.4/1.45 MM SHORT

JuggerKnot 1.4 mm Short with MaxBraid Suture

Part Number	Description	Suture Color	Needles	QTY
912068	Single Loaded #1 MB with needles	Blue/White	Yes	1
912069	Single Loaded #1 MB (10 pack)	Blue/White	Yes	10
912075	Implant System (Implant, Guide, Drill Bit)	Blue/White	Yes	

JuggerKnot 1.4 mm Short with BroadBand Tape

Part Number	Description	Suture Color	Needles	QTY
CM-99114BN	Single Loaded #1 (1.0 mm) BB	White/Blue	Yes	1
CM-99114BNT	Single Loaded #1 (1.0 mm) BB (10 pack)	White/Blue	Yes	10

JuggerKnot 1.4 mm Short Disposable Instrumentation

Part Number	Description
912071	1.4 mm Short Drill Bit
912073	1.4 mm Short Disposable Kit (Guide, Drill Bit)

JuggerKnot 1.4 mm Short Reusable Instrumentation

Part Number	Description
912072	1.4 mm Short Reusable Guide

JuggerKnot 1.45 mm Rigid with Maxbraid Suture


Part Number	Description	Suture Color	Needles	QTY
110005307	Single Loaded #1 MB (Includes Drill Bit)	Blue/White	Yes	1
110005319	Single Loaded #1 MB (10 pack)	Blue/White	Yes	10
110005315	Single Loaded #2 MB (Includes Drill Bit)	Blue	Yes	1
110005320	Single Loaded #2 MB (10 pack)	Blue	Yes	10

JuggerKnot 1.45 mm Rigid with BroadBand Tape

Part Number	Description	Suture Color	Needles	QTY
CM-99145BRN	Single Loaded #2 (1.5 mm) BB	Blue/Black	Yes	1

JuggerKnot 1.45 mm Rigid Disposable Instrumentation

Part Number	Description
110005323	1.45 mm Rigid Drill Bits (5 pack)




1.0 MM MINI

1.45 MM LONG

JUGGERKNOT 1.0 MM MINI AND 1.45 MM LONG ALL-SUTURE ANCHORS


The JuggerKnot 1.0 mm Mini Anchor is currently the smallest all-suture anchor on the market and was designed for hand and wrist procedures. It is loaded with either a 2-0 or 3-0 MaxBraid Suture with T-28 tapered needles.

The JuggerKnot 1.45 mm Long Soft Anchor is designed specifically for hip labral repair with an inserter 2.3 inches longer than the original 1.4 mm anchor to provide sufficient amount of working length in all areas of the hip labrum. It is offered single loaded with either a #1 or #2 MaxBraid Suture or a #2 (1.5 mm wide) BroadBand Tape.



1.0 MM MINI

**1.0 mm Mini with
MaxBraid Suture**



**1.45 mm Long with
MaxBraid Suture**



**1.45 mm Long with
BroadBand Tape**



JUGGERKNOT 1.0/1.45 mm LONG

JuggerKnot 1.0 mm Mini with Maxbraid Suture

Part Number	Description	Suture Color	Needles	QTY
912076	Single Loaded 2-0 MB (Includes Drill Bit)	White/Blue	Yes	1
912080	Single Loaded 2-0 MB (10 pack)	White/Blue	Yes	10
912082	Single Loaded 3-0 MB (Includes Drill Bit)	White/Blue	Yes	1
912084	Single Loaded 3-0 MB (10 pack)	White/Blue	Yes	10

JuggerKnot 1.0 mm Mini Disposable Instrumentation

Part Number	Description	Suture Color	Needles	QTY
912077	1.0 mm Mini Step Drill Bit	-	-	-

JuggerKnot 1.45 mm Long with Maxbraid Suture

Part Number	Description	Suture Color	Needles	QTY
916030	Single Loaded #1 MB	White	-	1
010000334	Single Loaded #1 MB (10 pack)	White	-	10
010000324	Single Loaded #2 MB	Blue	-	1
010000335	Single Loaded #2 MB (10 pack)	Blue	-	10

JuggerKnot 1.45 mm Long with BroadBand Tape


Part Number	Description	Suture Color	Needles	QTY
CM-991145BL	Single Loaded #2 (1.5 mm) BB	White/Blue	-	1
CM-992145BL	Single Loaded #2 (1.5 mm) BB	Black/White	-	1


JuggerKnot 1.45 mm Long Disposable Instrumentation

Part Number	Description
110016992	1.45 mm Long Flex Drill Bit with Centering Sleeve
916020	1.45 mm Long Curved Kit (Crvd Guide, Cntrng Slv, Flex Drill)
916040	1.45 mm Long Straight Kit (Straight Guide, Drill Bit)
916036	1.45 mm Long Flex Drill Bit
916037	1.45 mm Long Rigid Drill Bit

JuggerKnot 1.45 mm Long Reusable Instrumentation

Part Number	Description
110002759	1.45 mm Long Fishmouth Guide
916034	1.45 mm Long Curved Guide
916035	1.45 mm Long Straight Pronged Guide
916038	1.45 mm Long Trocar
916039	1.45 mm Long Obturator





JUGGERKNOTLESS 2.1 MM ALL-SUTURE ANCHOR

The 2.1 mm JuggerKnotless was the first all-suture knotless anchor on the market and was designed specifically for shoulder labral repair. It combines the innovation and clinical experience of the JuggerKnot All-Suture Anchor platform and the ZipLoop™ Technology.

The design of the implant eliminates the need for a surgeon to tie arthroscopic surgical knots. The locking mechanism is created in the design and technique itself and is pushed away from the articulating surface.

2.1 JuggerKnotless with MaxBraid Suture



JUGGERKNOTLESS 2.1 MM

JuggerKnotless 2.1 mm with MaxBraid Suture


Part Number	Description	Suture Color	QTY
110005198	Single Loaded	-	1
110005199	Single Loaded (10 pack)	-	10

JuggerKnotless 2.1 mm Disposable Instrumentation

Part Number	Description
110003172	JuggerKnotless Straight Kit (Guide, Drill Bit)
110003182	2.1 mm Drill Bit

JuggerKnotless 2.1 mm Reusable Instrumentation

Part Number	Description
110003181	JuggerKnotless Four-Prong Guide
110017436	JuggerKnotless Fishmouth Guide
110017434	JuggerKnotless Obturator
110017433	JuggerKnotless Trocar



SURELOCK 1.4 MM AND 2.2 MM ALL-SUTURE ANCHORS

The SureLock All-Suture Anchor family was designed to achieve consistent anchor delivery with an inserter-controlled deployment method delivering the anchor below the bone cortex. The SureLock Anchor family provides versatility in this comprehensive all-suture anchor portfolio and allows for confidence in anchor fixation.

The 1.4 mm SureLock Anchor is single loaded with #2 MaxBraid Suture, while the 2.2 mm SureLock Anchor is double loaded with #2 MaxBraid Suture.


**SureLock 1.4 mm with
#2 MaxBraid Suture**



**SureLock 2.2 mm with
#2 MaxBraid Suture**



**MECHANICAL TENSIONING
DESIGNED TO CONTROL
ANCHOR DEPLOYMENT
BELOW THE BONE CORTEX**



SURELOCK ANCHORS

SureLock 1.4 mm with MaxBraid Suture

Part Number	Description	Suture Color	QTY
CM-9614F	Single Loaded #2 MB	Blue/White	1

SureLock 1.4 mm Disposable Instrumentation

Part Number	Description
CM-9640	1.4 mm SureLock Drill Bit
CM-9640S	1.4 mm SureLock Drill Bit with Centering Sleeve

SureLock 1.4 mm Reusable Instrumentation

Part Number	Description
CM-9643	1.4 mm SureLock Straight Guide
CM-9645	1.4 mm SureLock Curved Guide
CM-9641	1.4 mm SureLock Obturator

SureLock 2.2 mm with MaxBraid Suture

Part Number	Description	Suture Color	QTY
CM-9622F	Double Loaded #2 MB	Blue/White, Black/White	1

SureLock 2.2 mm Disposable Instrumentation


Part Number	Description
CM-9620	2.2 mm SureLock Drill Bit
CM-9620S	2.2 mm SureLock Drill Bit with Centering Sleeve


SureLock 2.2 mm Reusable Instrumentation

Part Number	Description
CM-9623	2.2 mm SureLock Straight Guide
CM-9625	2.2 mm SureLock Curved Guide
CM-9621	2.2 mm SureLock Obturator



THE
LEADER
IN
ALL-SUTURE
ANCHOR
INNOVATION







3587.1-GLBL-en-Issue Date 2021-10-01

All content herein is protected by copyright, trademarks and other intellectual property rights, as applicable, owned by or licensed to Zimmer Biomet or its affiliates unless otherwise indicated, and must not be redistributed, duplicated or disclosed, in whole or in part, without the express written consent of Zimmer Biomet.

This material is intended for health care professionals. Distribution to any other recipient is prohibited.

For indications, contraindications, warnings, precautions, potential adverse effects and patient counseling information, see the package insert or contact your local representative; visit www.zimmerbiomet.com for additional product information.

Check for country product clearances and reference product specific instructions for use.

Not for distribution in France.

©2021 Zimmer Biomet



**JuggerKnot Soft Anchors
Legal Manufacturer:**

Biomet Orthopedics
56 East Bell Drive
P.O. Box 587
Warsaw, Indiana 46581 USA



**JuggerKnot with BroadBand
Legal Manufacturer:**

Riverpoint Medical
825 NE 25th Avenue,
Portland, OR 97232 USA

www.zimmerbiomet.com
www.rpmed.com