

# SpeedSnare™ Suture Passer System


First single step system with right and left bends.

## Not too big, not too small

- Effective combination of sturdiness and size
- 1.8 mm shaft to resist breakage
- Minimally invasive compared to standard 3.0 mm bird beak instrumentation

## Some like it locked, some don't

- The first system designed with “suture locked” and “suture sliding” positions
- Pass and retrieve in one step (accepts suture and suture tape)
- Audible and tactile feedback for confidence in snare positions



ZIMMER BIOMET  
Your progress. Our promise.®



## Ordering Information

Description	Type	Part Number
SpeedSnare Suture Passer	30° Up	110030129
	60° Up	110030130
	30° Right	110030131
	30° Left	110030132
	Straight	110030133

All content herein is protected by copyright, trademarks and other intellectual property rights, as applicable, owned by or licensed to Zimmer Biomet or its affiliates unless otherwise indicated, and must not be redistributed, duplicated or disclosed, in whole or in part, without the express written consent of Zimmer Biomet.

This material is intended for health care professionals. Distribution to any other recipient is prohibited.

For product information, including indications, contraindications, warnings, precautions, potential adverse effects and patient counseling information, see the package insert and [www.zimmerbiomet.com](http://www.zimmerbiomet.com).

Check for country product clearances and reference product specific instructions for use.

© 2017 Zimmer Biomet



**ZIMMER BIOMET**  
Your progress. Our promise.®



**Legal Manufacturer**  
T.A.G Medical Products  
Kibbutz Gaaton 25130  
Israel

[zimmerbiomet.com](http://zimmerbiomet.com)

**Distributor**

Biomet Sports Medicine  
56 East Bell Drive  
P.O. Box 587  
Warsaw, Indiana 46581  
USA